

# Lab Settings

24/10/2024 12:27 pm BST

## Summary

Lab Settings can be accessed via: Administration > Settings and selecting the 'Lab' category.

This page documents each setting.

**Configuration Settings [User - jupiter : Site - Esher : Profile - Back Reception]**

**Categories**

- Account
- Reporting
- Diary
- Reception
- Treatment
- Referral
- Lab**
- Vet Scan Lab
- Labelling
- SMS
- Email
- DayLists
- Svstem
- Miscellaneous
- HTML Labellir
- Smart Flow
- VetCheck
- Follow Up
- GDPR
- PACS
- Message

Site:

**Lab**

**External Lab**

**IDEXX Interlink Communication**

Service Host:

Service Port:

Client Port:

Client File Port:

File Download Timeout(seconds):

Connection Type:

Server URL:

Filtered Application Names (CSV):

For Idexx InterLink in Windows Terminal Session, please set INDIRECT TCP Connection Type and correct Server URL, example: http://[ServerNameOrIP][:PortOnlyIfRequired]/voyager-sessions-ejb/LabService?wsdl

**HESKA Data Capture Utility (DCU)**

DCU Output Folder Location:

Please make sure the folder is shared (Read/Write Access) for all Merlin Users.

☒ Filter By Client Category

## Idexx Interlink Communication.

These fields are used by MWI/Idexx to configure Idexx Interlink.

## HESKA Data Capture Utility (DCU)

As above, these fields are used by MWI/Idexx to configure Idexx Interlink.

## Filter by Client Category

**Enabled** = Allows users to filter Lab Results by Client Category in Utilities > General Lab.

**Disabled** = Do not allow users to filter Lab Results by Client Category in Utilities > General Lab.

