

Vet Scan Lab Settings


24/10/2024 12:29 pm BST

Summary

Vet Scan Settings can be accessed via: Administration > Settings and selecting the 'Vet Scan Lab' category.

This page is used to configure Abaxis Vet Scan Lab by our Engineers. This page documents each setting.

Details



The screenshot shows a web application interface titled "Configuration Settings [User - jupiter : Site - Esher : Profile - Back Reception]". On the left is a "Categories" sidebar with icons and labels for various settings: Account, Reporting, Diary, Reception, Treatment, Referral, Lab, Vet Scan Lab (highlighted), Labelling, SMS, Email, DayLists, System, Miscellaneous, HTML Labelling, Smart Flow, VetCheck, Follow Up, GDPR, PACS, and Message. The main content area is for the "Vet Scan Lab" category. At the top, there's a "Site" dropdown menu set to "Esher" and a "Save" button. Below this are four input fields: "Vet Scan URL" with the value "https://fuse03.abaxistest.com:8443/vetsync/v1", "Vet Scan User" with the value "pms", "Vet Scan Password" with masked characters, and "Check results every" with a value of "0" and the unit "seconds". At the bottom of the main area are two buttons: "Test Connection" and "Disable Result Checking (0 seconds)".

Vet Scan URL

The URL of the Vet Scan Fuse box.

Vet Scan User

Enter username here.

Vet Scan Password

Enter password here.

Check results every X seconds

Controls how often the system checks for results.

